

PKLR antibody

Product Information

Catalog No.:	FNab06490
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a key role in glycolysis.

Immunogen information


Immunogen:	pyruvate kinase, liver and RBC
Synonyms:	Pyruvate kinase PKLR Pyruvate kinase 1 Pyruvate kinase isozymes L/R R-type/L-type pyruvate kinase Red cell/liver pyruvate kinase PKLR PK1 PKL
Observed MW:	58-62 kDa
Uniprot ID :	P30613

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab06490(PKLR Antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab06490(PKLR antibody) at dilution of 1:1000