

PKNOX2 antibody

Product Information

Catalog No.:	FNab06500
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	PBX/knotted 1 homeobox 2
Synonyms:	Homeobox protein PKNOX2 Homeobox protein PREP-2 PBX/knotted homeobox 2 PKNOX2 PREP2
Observed MW:	70 kDa, 55 kDa
Uniprot ID :	Q96KN3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab06500(PKNOX2 antibody) at dilution of 1:50



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab06500(PKNOX2 antibody) at dilution of 1:300