

PLA2G4F antibody

Product Information

Catalog No.:	FNab06507
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Calcium-dependent phospholipase A2 that selectively hydrolyzes glycerophospholipids in the sn-2 position. Has higher enzyme activity for phosphatidylethanolamine than phosphatidylcholine(By similarity).

Immunogen information

Immunogen:	phospholipase A2, group IVF
Synonyms:	Cytosolic phospholipase A2 zeta (cPLA2-zeta) Phospholipase A2 group IVF PLA2G4F
Observed MW:	None
Uniprot ID :	Q68DD2

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNab06507(PLA2G4F antibody) at dilution of 1:200