

---

## PLAA antibody

### Product Information

Catalog No.:	FNab06508
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Involved in the maintenance of ubiquitin levels.

### Immunogen information


Immunogen:	phospholipase A2-activating protein
Synonyms:	Phospholipase A-2-activating protein (PLA2P, PLAP) PLAA PLAP
Observed MW:	72-74 kDa
Uniprot ID :	Q9Y263

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-PLAA (IP:FNab06508, 4ug;  
Detection:FNab06508 1:600) with K-562 cells  
lysate 3200ug.



K-562 cells were subjected to SDS PAGE  
followed by western blot with FNab06508(PLAA  
antibody) at dilution of 1:500