

## PLAGL1 antibody

### Product Information

Catalog No.:	FNab06516
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Shows weak transcriptional activatory activity. Transcriptional regulator of the type 1 receptor for pituitary adenylate cyclase-activating polypeptide.

### Immunogen information

Immunogen:	pleiomorphic adenoma gene-like 1
Synonyms:	Zinc finger protein PLAGL1 Lost on transformation 1 (LOT-1) Pleiomorphic adenoma-like protein 1 Tumor suppressor ZAC PLAGL1 LOT1 ZAC
Observed MW:	51 kDa
Uniprot ID :	Q9UM63

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab06516( PLAGL1 Antibody) at dilution of 1:600