

PLLP antibody

Product Information

Catalog No.:	FNab06524
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	plasma membrane proteolipid(plasmolipin)
Synonyms:	Plasmolipin Plasma membrane proteolipid PLLP PMLP TM4SF11
Observed MW:	18 kDa
Uniprot ID :	Q9Y342

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:1000-1:4000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab06524(Plasmolipin antibody) at dilution of 1:1000