

PLEK2 antibody

Product Information

Catalog No.:	FNab06533
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May help orchestrate cytoskeletal arrangement. Contribute to lamellipodia formation.

Immunogen information


Immunogen:	pleckstrin 2
Synonyms:	Pleckstrin-2 PLEK2
Observed MW:	38-40 kDa
Uniprot ID :	Q9NYT0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab06533(PLEK2 antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab06533(PLEK2 antibody) at dilution of 1:500