

PLEKHA3 antibody

Product Information

Catalog No.:	FNab06535
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in Golgi to cell surface membrane traffic. Induces membrane tubulation. Binds preferentially to phosphatidylinositol 4-phosphate(PtdIns4P).

Immunogen information

Immunogen:	pleckstrin homology domain containing, family A(phosphoinositide binding specific) member 3
Synonyms:	Pleckstrin homology domain-containing family A member 3 (PH domain-containing family A member 3) Phosphatidylinositol-four-phosphate adapter protein 1 (FAPP-1, Phosphoinositol 4-phosphate adapter protein 1) PLEKHA3 FAPP1
Observed MW:	34 kDa
Uniprot ID :	Q9HB20

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab06535 (PLEKHA3 antibody) at dilution of 1:200