

PLEKHA5 antibody

Product Information

Catalog No.:	FNab06536
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	pleckstrin homology domain containing, family A member 5
Synonyms:	Pleckstrin homology domain-containing family A member 5 (PH domain-containing family A member 5) Phosphoinositol 3-phosphate-binding protein 2 (PEPP-2) PLEKHA5 KIAA1686 PEPP2
Observed MW:	127 kDa, 74 kDa
Uniprot ID :	Q9HAU0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab06536(PLEKHA5 antibody) at dilution of 1:500