

PLEKHH2 antibody

Product Information

Catalog No.:	FNab06542
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


In the kidney glomerulus may play a role in linking podocyte foot processes to the glomerular basement membrane. May be involved in stabilization of F-actin by attenuating its depolymerization. Can recruit TGFB1I1 from focal adhesions to podocyte lamellipodia.

Immunogen information

Immunogen:	pleckstrin homology domain containing, family H(with MyTH4 domain) member 2
Synonyms:	Pleckstrin homology domain-containing family H member 2 PLEKHH2 KIAA2028
Observed MW:	120 kDa
Uniprot ID :	Q8IVE3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab06542(PLEKHH2 Antibody) at dilution of 1:300