

PLEKHO2 antibody

Product Information

Catalog No.:	FNab06544
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	pleckstrin homology domain containing, family O member 2
Synonyms:	Pleckstrin homology domain-containing family O member 2 (PH domain-containing family O member 2) Pleckstrin homology domain-containing family Q member 1 (PH domain-containing family Q member 1) PLEKHO2 PLEKHQ1
Observed MW:	53 kDa
Uniprot ID :	Q8TD55

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung tissue slide using FNab06544(PLEKHO2 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06544(PLEKHO2 Antibody) at dilution of 1:300