

## PLS3 antibody

### Product Information

Catalog No.:	FNab06559
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Actin-bundling protein found in intestinal microvilli, hair cell stereocilia, and fibroblast filopodia. May play a role in the regulation of bone development.

### Immunogen information


Immunogen:	plastin 3(T isoform)
Synonyms:	Plastin-3 T-plastin PLS3
Observed MW:	63 kDa, 38 kDa
Uniprot ID :	P13797

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using FNab06559(PLS3 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06559(PLS3 antibody) at dilution of 1:500