

## PLSCR4 antibody

### Product Information

Catalog No.:	FNab06562
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May mediate accelerated ATP-independent bidirectional transbilayer migration of phospholipids upon binding calcium ions that results in a loss of phospholipid asymmetry in the plasma membrane. May play a central role in the initiation of fibrin clot formation, in the activation of mast cells and in the recognition of apoptotic and injured cells by the reticuloendothelial system.

### Immunogen information

Immunogen:	phospholipid scramblase 4
Synonyms:	Phospholipid scramblase 4 (PL scramblase 4) Ca(2+)-dependent phospholipid scramblase 4 Cell growth-inhibiting gene 43 protein TRA1 PLSCR4
Observed MW:	37 kDa
Uniprot ID :	Q9NRQ2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human heart tissue were subjected to SDS PAGE followed by western blot with FNab06562(PLSCR4 antibody) at dilution of 1:400