

PLVAP antibody

Product Information

Catalog No.:	FNab06564
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the formation of stomatal and fenestral diaphragms of caveolae. May function in microvascular permeability.

Immunogen information

Immunogen:	plasmalemma vesicle associated protein
Synonyms:	Plasmalemma vesicle-associated protein Fenestrated endothelial-linked structure protein Plasmalemma vesicle protein 1 (PV-1) PLVAP FELS PV1
Observed MW:	51 kDa
Uniprot ID :	Q9BX97

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB : 1:500-1:1000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06564(PLVAP antibody) at dilution of 1:800