

## ATP2B2 antibody

### Product Information

Catalog No.:	FNab06568
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	ATPase, Ca <sup>++</sup> transporting, plasma membrane 2
Synonyms:	Plasma membrane calcium-transporting ATPase 2 (PMCA2) Plasma membrane calcium ATPase isoform 2 Plasma membrane calcium pump isoform 2 ATP2B2 PMCA2
Observed MW:	137 kDa
Uniprot ID :	Q01814

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06568(ATP2B2 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab06568(ATP2B2 antibody) at dilution of 1:500