

PMFBP1 antibody

Product Information

Catalog No.:	FNab06571
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in sperm morphology especially the sperm tail and consequently affect fertility. May also be involved in the general organization of cellular cytoskeleton.

Immunogen information

Immunogen:	polyamine modulated factor 1 binding protein 1
Synonyms:	Polyamine-modulated factor 1-binding protein 1 (PMF-1-binding protein) PMFBP1
Observed MW:	119 kDa
Uniprot ID :	Q8TBY8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab06571(PMFBP1 antibody) at dilution of 1:500