

PNMA2 antibody

Product Information

Catalog No.:	FNab06586
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	paraneoplastic antigen MA2
Synonyms:	Paraneoplastic antigen Ma2 40 kDa neuronal protein Onconeuroal antigen Ma2 Paraneoplastic neuronal antigen MM2 PNMA2 KIAA0883 MA2
Observed MW:	40 kDa
Uniprot ID :	Q9UL42

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab06586(PNMA2 antibody) at dilution of 1:2000