

PNMA8A antibody

Product Information

Catalog No.:	FNab06587
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	PNMA-like 1
Synonyms:	Paraneoplastic antigen-like protein 8A PNMA-like protein 1 PNMA8A PNMAL1
Observed MW:	48 kDa
Uniprot ID :	Q86V59

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06587(PNMAL1 antibody) at dilution of 1:500