

## PODXL2 antibody

### Product Information

Catalog No.:	FNab06600
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Acts as a ligand for vascular selectins. Mediates rapid rolling of leukocytes over vascular surfaces through high affinity divalent cation-dependent interactions with E-, P-and L-selectins.

### Immunogen information

Immunogen:	podocalyxin-like 2
Synonyms:	Podocalyxin-like protein 2 Endoglycan PODXL2
Observed MW:	65 kDa
Uniprot ID :	Q9NZ53

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab06600(PODXL2 antibody) at dilution of 1:400