

POF1B antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab06601 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Plays a key role in the organization of epithelial monolayers by regulating the actin cytoskeleton. May be involved in ovary development.

Immunogen information

| | |
|--------------|--|
| Immunogen: | premature ovarian failure, 1B |
| Synonyms: | Protein POF1B Premature ovarian failure protein 1B POF1B |
| Observed MW: | 69 kDa |
| Uniprot ID : | Q8WVV4 |

Application

| | |
|-----------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution: | WB: 1:500-1:2000 |
| Image: | |



human bladder tissue were subjected to SDS
PAGE followed by western blot with
FNab06601(POF1B antibody) at dilution of 1:500