

## POLDIP2 antibody

### Product Information

Catalog No.:	FNab06609
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	polymerase(DNA-directed), delta interacting protein 2
Synonyms:	Polymerase delta-interacting protein 2 38 kDa DNA polymerase delta interaction protein (p38) POLDIP2 PDIP38 POLD4
Observed MW:	38 kDa
Uniprot ID :	Q9Y2S7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06609(POLDIP2 antibody) at dilution of 1:500