

## POTEA antibody

### Product Information

Catalog No.:	FNab06646
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	POTE ankyrin domain family, member A
Synonyms:	POTE ankyrin domain family member A ANKRD26-like family A member 1 Prostate, ovary, testis-expressed protein on chromosome 8 (POTE-8) POTEA A26A1 POTE8
Observed MW:	53 kDa
Uniprot ID :	Q6S8J7

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000 ; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab06646(POTEA Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab06646(POTEA Antibody) at dilution of 1:300