

POTEA-Specific antibody

Product Information

Catalog No.:	FNab06647
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	POTE ankyrin domain family, member A
Synonyms:	POTE ankyrin domain family member A ANKRD26-like family A member 1 Prostate, ovary, testis-expressed protein on chromosome 8 (POTE-8) POTEA A26A1 POTE8
Observed MW:	None
Uniprot ID :	Q6S8J7

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, IF
Recommended dilution:	IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemical of paraffin-embedded human ovary using FNab06647(POTEA-Specific antibody) at dilution of 1:50