

PTPA antibody

Product Information

Catalog No.:	FNab06655
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	protein phosphatase 2A activator, regulatory subunit 4
Synonyms:	Serine/threonine-protein phosphatase 2A activator PP2A, subunit B', PR53 isoform Phosphotyrosyl phosphatase activator (PTPA) Serine/threonine-protein phosphatase 2A regulatory subunit 4 Serine/threonine-protein phosphatase 2A regulatory subunit B' PTPA PPP2R4
Observed MW:	37 kDa
Uniprot ID :	Q15257

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000 ; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab06655(PP2A Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06655(PP2A antibody) at dilution of 1:500