

PLPP7 antibody

Product Information

Catalog No.:	FNab06659
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	phosphatidic acid phosphatase type 2 domain containing 3
Synonyms:	Inactive phospholipid phosphatase 7 Phosphatidic acid phosphatase type 2 domain-containing protein 3 PLPP7 C9orf67 PPAPDC3
Observed MW:	29 kDa
Uniprot ID :	Q8NBV4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab06659(PPAPDC3 antibody) at dilution of 1:800