

PPARG antibody

Product Information

Catalog No.:	FNab06661
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	peroxisome proliferator-activated receptor gamma
Synonyms:	Peroxisome proliferator-activated receptor gamma (PPAR-gamma) Nuclear receptor subfamily 1 group C member 3 PPARG NR1C3
Observed MW:	50-60 kDa
Uniprot ID :	P37231

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab06661(PPARG Antibody) at dilution of 1:200



K-562 cells were subjected to SDS PAGE followed by western blot with FNab06661(PPARG antibody) at dilution of 1:1000