

PPM1L antibody

Product Information

Catalog No.:	FNab06689
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as a suppressor of the SAPK signaling pathways by associating with and dephosphorylating MAP3K7/TAK1 and MAP3K5, and by attenuating the association between MAP3K7/TAK1 and MAP2K4 or MAP2K6.

Immunogen information


Immunogen:	protein phosphatase 1(formerly 2C)-like
Synonyms:	Protein phosphatase 1L Protein phosphatase 1-like Protein phosphatase 2C isoform epsilon (PP2C-epsilon) PPM1L PP2CE
Observed MW:	41-45 kDa
Uniprot ID :	Q5SGD2

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab06689(PPM1L Antibody) at dilution of 1:50



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab06689(PPM1L antibody) at dilution of 1:500