

PPME1 antibody

Product Information

Catalog No.:	FNab06692
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Demethylates proteins that have been reversibly carboxymethylated. Demethylates PPP2CB(in vitro) and PPP2CA. Binding to PPP2CA displaces the manganese ion and inactivates the enzyme.

Immunogen information

Immunogen:	protein phosphatase methylesterase 1
Synonyms:	Protein phosphatase methylesterase 1 (PME-1) PPME1 PME1
Observed MW:	44 kDa
Uniprot ID :	Q9Y570

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06692(PPME1 antibody) at dilution of 1:1000