

PPP1R12B antibody

Product Information

Catalog No.:	FNab06700
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Regulates myosin phosphatase activity. Augments Ca(2+) sensitivity of the contractile apparatus.

Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 12B
Synonyms:	Protein phosphatase 1 regulatory subunit 12B Myosin phosphatase-targeting subunit 2 (Myosin phosphatase target subunit 2) PPP1R12B MYPT2
Observed MW:	110 kDa
Uniprot ID :	O60237

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab06700(PPP1R12B antibody) at dilution of 1:10000