

PPP1R14B antibody

Product Information

Catalog No.:	FNab06703
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibitor of PPP1CA. Has over 50-fold higher inhibitory activity when phosphorylated(By similarity).

Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 14B
Synonyms:	Protein phosphatase 1 regulatory subunit 14B Phospholipase C-beta-3 neighbouring gene protein PPP1R14B PLCB3N PNG
Observed MW:	23 kDa
Uniprot ID :	Q96C90

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab06703(PPP1R14B antibody) at dilution of 1:300