

## PPP1R7 antibody

### Product Information

Catalog No.:	FNab06708
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Regulatory subunit of protein phosphatase 1.

### Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 7
Synonyms:	Protein phosphatase 1 regulatory subunit 7 Protein phosphatase 1 regulatory subunit 22 PPP1R7 SDS22
Observed MW:	42 kDa
Uniprot ID :	Q15435

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab06708(PPP1R7 antibody) at dilution of 1:1000