

PPP3CC antibody

Product Information

Catalog No.:	FNab06729
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Calcium-dependent, calmodulin-stimulated protein phosphatase. This subunit may have a role in the calmodulin activation of calcineurin.

Immunogen information

Immunogen:	protein phosphatase 3(formerly 2B), catalytic subunit, gamma isoform
Synonyms:	Serine/threonine-protein phosphatase 2B catalytic subunit gamma isoform CAM-PRP catalytic subunit Calcineurin, testis-specific catalytic subunit Calmodulin-dependent calcineurin A subunit gamma isoform PPP3CC CALNA3 CNA3
Observed MW:	57 kDa
Uniprot ID :	P48454

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab06729(PPP3CC antibody) at dilution of 1:500