

DESI2 antibody

Product Information

Catalog No.:	FNab06737
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	PPPDE peptidase domain containing 1
Synonyms:	Deubiquitinase DESI2 Desumoylating isopeptidase 2 (DeSI-2) PPPDE peptidase domain-containing protein 1 Protein FAM152A DESI2 C1orf121 FAM152A PPPDE1
Observed MW:	23 kDa
Uniprot ID :	Q9BSY9

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer using FNaB06737(PPPDE1 antibody) at dilution of 1:50



IP Result of anti-PPPDE1/PNAS4 (IP:FNaB06737, 4ug; Detection:FNaB06737 1:500) with A549 cells lysate 2800ug.



A549 cells were subjected to SDS PAGE followed by western blot with FNaB06737(PPPDE1 antibody) at dilution of 1:500