

## PPRC1 antibody

### Product Information

Catalog No.:	FNab06738
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Acts as a coactivator during transcriptional activation of nuclear genes related to mitochondrial biogenesis and cell growth. Involved in the transcription coactivation of CREB and NRF1 target genes.

### Immunogen information

Immunogen:	peroxisome proliferator-activated receptor gamma, coactivator-related 1
Synonyms:	Peroxisome proliferator-activated receptor gamma coactivator-related protein 1 PGC-1-related coactivator (PRC) PPRC1 KIAA0595
Observed MW:	177 kDa
Uniprot ID :	Q5VV67

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200 - 1:2000
Image:	



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab06738( PPRC1 Antibody) at dilution of 1:300