

## PPT1 antibody

### Product Information

Catalog No.:	FNab06739
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Removes thioester-linked fatty acyl groups such as palmitate from modified cysteine residues in proteins or peptides during lysosomal degradation. Prefers acyl chain lengths of 14 to 18 carbons.

### Immunogen information

Immunogen:	palmitoyl-protein thioesterase 1
Synonyms:	Palmitoyl-protein thioesterase 1 (PPT-1) Palmitoyl-protein hydrolase 1 PPT1 CLN1 PPT
Observed MW:	None
Uniprot ID :	P50897

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab06739(PPT1 antibody) at dilution of 1:50