

## PPWD1 antibody

### Product Information

Catalog No.:	FNab06741
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Putative peptidylprolyl isomerase(PPIase). PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. May be involved in pre-mRNA splicing.

### Immunogen information

Immunogen:	peptidylprolyl isomerase domain and WD repeat containing 1
Synonyms:	Peptidylprolyl isomerase domain and WD repeat-containing protein 1 Spliceosome-associated cyclophilin PPWD1 KIAA0073
Observed MW:	65-66 kDa
Uniprot ID :	Q96BP3

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



Sp2/0 cells were subjected to SDS PAGE followed by western blot with FNab06741(PPWD1 antibody) at dilution of 1:1000