

PRDM15 antibody

Product Information

Catalog No.:	FNab06753
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in transcriptional regulation.

Immunogen information


Immunogen:	PR domain containing 15
Synonyms:	PR domain zinc finger protein 15 PR domain-containing protein 15 Zinc finger protein 298 PRDM15 C21orf83 ZNF298
Observed MW:	169 kDa
Uniprot ID :	P57071

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:200-1:1000; IP: 1:500-1:1000
Image:	



IP Result of anti-PRDM15 (IP:FNab06753, 5ug;
Detection:FNab06753 1:500) with MCF-7 cells
lysate 2500ug.



MCF-7 cells were subjected to SDS PAGE
followed by western blot with
FNab06753(PRDM15 Antibody) at dilution of
1:300