

PRH1 antibody

Product Information

Catalog No.:	FNab06772
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	proline-rich protein HaeIII subfamily 1
Synonyms:	Salivary acidic proline-rich phosphoprotein 1/2 Db-s PRP-1/PRP-2 Parotid acidic protein (Pa) Parotid double-band protein Parotid isoelectric focusing variant protein (PIF-S) Parotid proline-rich protein 1/2 Pr1/Pr2 Protein C Salivary acidic proline-rich phosphoprotein 1/2 Salivary acidic proline-rich phosphoprotein 3/4 Alternative names: Db-F PIF-F PRP-3/PRP-4 Protein A Peptide P-C PRH1
Observed MW:	17 kDa
Uniprot ID :	P02810

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IF: 1:20-1:200
Image:	

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06772(PRH1 Antibody) at dilution of 1:300

