

## PRICKLE3 antibody

### Product Information

Catalog No.:	FNab06775
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	prickle homolog 3(Drosophila)
Synonyms:	Prickle planar cell polarity protein 3 LIM domain only protein 6 (LMO-6) Prickle-like protein 3 (Pk3) Triple LIM domain protein 6 PRICKLE3 LMO6
Observed MW:	69 kDa
Uniprot ID :	O43900

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab06775(PRICKLE3 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab06775(PRICKLE3 antibody) at dilution of 1:500