

PHB1 antibody

Product Information

Catalog No.:	FNab06800
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5F2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	prohibitin
Synonyms:	Prohibitin 1 PHB1 PHB
Observed MW:	30 kDa
Uniprot ID :	P35232

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab06800(PHB antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06800(PHB antibody) at dilution of 1:500