

PROS1 antibody

Product Information

Catalog No.:	FNab06803
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Anticoagulant plasma protein; it is a cofactor to activated protein C in the degradation of coagulation factors Va and VIIIa. It helps to prevent coagulation and stimulating fibrinolysis.

Immunogen information


Immunogen:	protein S(alpha)
Synonyms:	Vitamin K-dependent protein S PROS1 PROS
Observed MW:	75 kDa
Uniprot ID :	P07225

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab06803(PROS1 antibody) at dilution of 1:100



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab06803(PROS1 antibody) at dilution of 1:1000