

LZTS3 antibody

Product Information

Catalog No.:	FNab06804
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ProSAPiP1 protein
Synonyms:	Leucine zipper putative tumor suppressor 3 ProSAP-interacting protein 1 (ProSAPiP1) LZTS3 KIAA0552 PROSAPIP1
Observed MW:	75 kDa
Uniprot ID :	O60299

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06804(ProSAPiP1 Antibody) at dilution of 1:600