

PLPBP antibody

Product Information

Catalog No.:	FNab06805
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	proline synthetase co-transcribed homolog(bacterial)
Synonyms:	Pyridoxal phosphate homeostasis protein (PLP homeostasis protein) Proline synthase co-transcribed bacterial homolog protein Pyridoxal phosphate-binding protein PLPBP PROSC
Observed MW:	30 kDa
Uniprot ID :	O94903

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab06805(PROSC antibody) at dilution of 1:3000