

## SERPINA5 antibody

### Product Information

Catalog No.:	FNab06806
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	serpin peptidase inhibitor, clade A(alpha-1 antiproteinase, antitrypsin), member 5
Synonyms:	Plasma serine protease inhibitor Acrosomal serine protease inhibitor Plasminogen activator inhibitor 3 (PAI-3, PAI3) Protein C inhibitor (PCI) Serpine A5 SERPINA5 PCI PLANH3 PROCI
Observed MW:	46 kDa
Uniprot ID :	P05154

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab06806(SERPINA5 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06806(SERPINA5 antibody) at dilution of 1:1000