

SERPINA5 antibody

Product Information

Catalog No.:	FNab06807
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2H5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	serpin peptidase inhibitor, clade A(alpha-1 antiproteinase, antitrypsin), member 5
Synonyms:	Plasma serine protease inhibitor Acrosomal serine protease inhibitor Plasminogen activator inhibitor 3 (PAI-3, PAI3) Protein C inhibitor (PCI) Serpin A5 SERPINA5 PCI PLANH3 PROCI
Observed MW:	46 kDa
Uniprot ID :	P05154

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab06807(Protein C inhibitor antibody) at dilution of 1:1000