

PRPF19 antibody

Product Information

Catalog No.:	FNab06813
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).

Immunogen information


Immunogen:	PRP19/PSO4 pre-mRNA processing factor 19 homolog
Synonyms:	Pre-mRNA-processing factor 19 Nuclear matrix protein 200 PRP19/PSO4 homolog (hPso4) RING-type E3 ubiquitin transferase PRP19 Senescence evasion factor PRPF19 NMP200 PRP19 SNEV
Observed MW:	60 kDa
Uniprot ID :	Q9UMS4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:100 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat brain using FNab06813(PRPF19 Antibody) at dilution of 1:100



SW620 cells were subjected to SDS PAGE followed by western blot with FNab06813(PRPF19 antibody) at dilution of 1:1000