

PRPF3 antibody

Product Information

Catalog No.:	FNab06814
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

Immunogen information


Immunogen:	PRP3 pre-mRNA processing factor 3 homolog
Synonyms:	U4/U6 small nuclear ribonucleoprotein Prp3 Pre-mRNA-splicing factor 3 (hPrp3) U4/U6 snRNP 90 kDa protein PRPF3 HPRP3 PRP3
Observed MW:	90 kDa or 78 kDa
Uniprot ID :	O43395

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-PRPF3 (IP:FNab06814, 5ug;
Detection:FNab06814 1:500) with mouse liver
tissue lysate 4000ug.



Y79 cells were subjected to SDS PAGE followed
by western blot with FNab06814(PRPF3 antibody)
at dilution of 1:1000