

PRPF39 antibody

Product Information

Catalog No.:	FNab06816
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in pre-mRNA splicing.

Immunogen information


Immunogen:	PRP39 pre-mRNA processing factor 39 homolog
Synonyms:	Pre-mRNA-processing factor 39 PRP39 homolog PRPF39
Observed MW:	75-85 kDa
Uniprot ID :	Q86UA1

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab06816 (PRPF39 Antibody) at dilution of 1:100



K-562 cells were subjected to SDS PAGE followed by western blot with FNab06816 (PRPF39 Antibody) at dilution of 1:1000