

PRPF40A antibody

Product Information

Catalog No.:	FNab06818
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Binds to WASL/N-WASP and suppresses its translocation from the nucleus to the cytoplasm, thereby inhibiting its cytoplasmic function(By similarity). Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape and migration. May play a role in cytokinesis. May be involved in pre-mRNA splicing.


Immunogen information

Immunogen:	PRP40 pre-mRNA processing factor 40 homolog A
Synonyms:	Pre-mRNA-processing factor 40 homolog A Fas ligand-associated factor 1 Formin-binding protein 11 Formin-binding protein 3 Huntingtin yeast partner A Huntingtin-interacting protein 10 (HIP-10) Huntingtin-interacting protein A Renal carcinoma antigen NY-REN-6 PRPF40A FBP11 FLAF1 FNBP3 HIP10 HYPA
Observed MW:	130 kDa
Uniprot ID :	O75400

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:2000

Image:



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06818(PRPF40A Antibody) at dilution of 1:1000