

PRPF6 antibody

Product Information

Catalog No.:	FNab06820
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in pre-mRNA splicing as component of the U4/U6-U5 tri-snRNP complex, one of the building blocks of the spliceosome. Enhances dihydrotestosterone-induced transactivation activity of AR, as well as dexamethasone-induced transactivation activity of NR3C1, but does not affect estrogen-induced transactivation.

Immunogen information


Immunogen:	PRP6 pre-mRNA processing factor 6 homolog
Synonyms:	Pre-mRNA-processing factor 6 Androgen receptor N-terminal domain-transactivating protein 1 (ANT-1) PRP6 homolog U5 snRNP-associated 102 kDa protein (U5-102 kDa protein) PRPF6 C20orf14
Observed MW:	100 kDa
Uniprot ID :	O94906

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab06820 (PRPF6 Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06820 (PRPF6 Antibody) at dilution of 1:1000